



PRECISION
VISION, LTD.
CONTACT LENS LABORATORY
Quality and Service Now....
With an Eye on the Future

1490 S. Yampa Way
Aurora, CO 80017
303.743.0494
800.843.5367
FAX 303.632.8216

www.precisionvision.com

Vertex Adjustment Table

MINUS	MANIFEST	PLUS	MINUS	MANIFEST	PLUS
3.87	4.00	4.25	9.87	11.25	13.00
4.00	4.25	4.50	10.12	11.50	13.37
4.25	4.50	4.75	10.25	11.75	13.75
4.50	4.75	5.00	10.50	12.00	14.00
4.75	5.00	5.25	10.62	12.25	14.37
5.00	5.25	5.62	10.87	12.50	14.75
5.12	5.50	5.87	11.12	12.75	15.00
5.37	5.75	6.12	11.25	13.00	15.50
5.62	6.00	6.50	11.50	13.25	15.75
5.75	6.25	6.75	11.62	13.50	16.12
6.00	6.50	7.00	11.87	13.75	16.50
6.25	6.75	7.37	12.00	14.00	16.75
6.50	7.00	7.62	12.12	14.25	17.25
6.62	7.25	8.00	12.37	14.50	17.50
6.87	7.50	8.25	12.50	14.75	18.00
7.12	7.75	8.50	12.75	15.00	18.25
7.25	8.00	8.87	12.87	15.25	18.75
7.50	8.50	9.50	13.00	15.50	19.00
7.87	8.75	9.75	13.25	15.75	19.37
8.12	9.00	10.12	13.50	16.00	19.75
8.37	9.25	10.37	13.50	16.25	20.25
8.62	9.50	10.75	13.75	16.75	21.00
8.75	9.75	11.00	14.12	17.00	21.50
8.87	10.00	11.37	14.25	17.25	21.75
9.12	10.25	11.75	14.50	17.50	22.25
9.37	10.50	12.00	14.62	17.75	22.50
9.50	10.75	12.37	14.75	18.00	23.00
9.75	11.00	12.75	15.00	18.25	23.50